



SINTEF

Micro-ElectroMechanical Systems (MEMS) Fabry-Perot (FP) spectrometer for measurement of contamination deposition

Authors: Grégory Bouquet¹, Kari Anne Bakke¹, Jon Tschudi¹, Runar Dahl-Hansen¹, and Julien Eck²

¹ Smart Sensor Systems and Microsystems, SINTEF, Norway

² Cleanliness and Contamination Control Engineer, Materials' Physics & Chemistry Section (TEC-QEE), ESTEC - European Space Research and Technology Centre, Keplerlaan 1, PO Box 299, 2200AG Noordwijk, The Netherlands

Acknowledgements: Atle Rognmo³, Hilde Isachsen³, Ingeborg Nævra Prestholdt³, Dung Do Dang³ and Thor Bakke³

³ Tunable AS, Kabelgata 8B, 0580 Oslo

ISMSE-ICPMSE 2024,
7-11 October 2024
Saint-Raphaël, France



SINTEF

Scope

- Background
- Standard measurement
- MEMS spectrometer development
- Some measurement results
- Conclusion



SINTEF

Background I

- Current monitoring methods of molecular contamination are based on witness exposure, allowing only the punctual and ex-situ estimation of the deposited levels on sensitive surfaces.
- Besides, conventional monitoring is performed usually until encapsulation, in the best case under fairing prior launch, and from launch to end of-life, only estimation by modelling is performed.
- Overall, there is currently a lack of representativeness and only a partial coverage of the molecular contamination monitoring over the whole mission.
- Having the capability to perform in-situ and real-time monitoring of the molecular contamination will improve the representativeness of the measurements and the triggering of mitigation actions to ensure optimal system performances until end-of-life (both on-ground and in-flight).



SINTEF

In situ monitoring of contamination

- In collaboration with European Space Agency, Tunable AS and SINTEF are aiming at developing a sensor breadboard at TRL 3 able to quantify the deposited molecular contamination (hydrocarbon, ester, and, silicones) to provide reliable in-situ and real-time data during assembly integration and testing with the possibility to be further developed for in-flight.
- The sensor is based on a piezo-MEMS Fabry-Perot Interferometer (FPI) developed and commercialized by Tunable AS for gas detection.
- The measurement system has been developed, mounted and tests on some highly concentrated reference samples initiated





SINTEF


Standard measurement of contamination⁽¹⁾



Extraction of peak height at a few characteristic wavelengths



Quantification using Beer Lambert law







⁽¹⁾ Detection of organic contamination of surfaces by infrared spectroscopy, ECSS-Q-ST-70-05C Rev.1 15 October 2019

⁽²⁾ Picture: courtesy of Lucideon <https://cdn.lucideon.com/v2/downloads/MOC-Technical-Poster.pdf>



SINTEF

Fabry-Perot spectrometer





a) FP filter; b) Lens L2; c) lens L1; d) Detector and detector front end electronic; e) Light source



- Volume: 15 cm x 15 cm x 10 cm.
- FWHM ranging between 25 cm^{-1} and 40 cm^{-1} (3)

(3) FWHM \propto path difference between in and out beam


Piezo MEMS Tunable filter



interference on thin oil film on water



micromirrors on wafer Photo: Karolina Milenko / SINTEF




Picture: courtesy of Tunable AS




SINTEF

FTIR versus Micro MEMS spectrometer




- A MEMS spectrometer would outperform FTIR in terms of volume, weight and cost, but with a trade off in term of spectral resolution and Signal to Noise Ratio =>no free lunch





- Using more points from a low-resolution spectra will enable good enough quantification, by using machine learning algorithm (e.g Partial Least Square) to handle collinearity between spectral channels



MEMS manufacturing and reliability



Current batch







SINTEF

A few examples of measurement

- Measurement performed using the MEMS spectrometer on two references samples
- Compared with FTIR measurement performed at 32 cm^{-1}
- Very good agreement between the two measurements showing that the FWHM is as expected



poly (dimethylsiloxane)



bis (2-ethylhexyl) phthalate



SINTEF

Conclusion

- In a project with ESA, we are aiming at developing a MEMS spectrometer for in-situ monitoring surface contamination for clean room and in-flight applications.
- We have developed a MEMS spectrometer using a MEMS Fabry Perot tunable filter at TRL 3 and performed tests showing that the spectrometer perform as required.
- The next step is to develop a stand alone unit allowing to collect contamination spectra in situ and establish calibration models for the four mains contaminants



SINTEF

Acknowledgement



Activity undertaken under ESA contract No. 4000139400

The Tunable logo, which is the word 'tunable' in a white, lowercase, sans-serif font, centered within a dark green rectangular background.

tunable

<https://www.tunable.com/>: Design the chips and developed and deposited the coatings, developed methods and equipment for assembling the chips, the electronics, calibrated the filters and developed the control algorithms.



SINTEF

Technology for a better society